## BT20 Biological Indicator

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**Quality certification**  
**Certificado de calidad**  
**Bionova® BT20**  
 Steam sterilization / Esterilización por Vapor  
*Geobacillus stearothermophilus* ATCC® 7953

**LOT**

Population / Población _____	CFU / UFC
D - value / Valor D (121 °C) _____	min.
Survival time / Tiempo de sobrevida _____	min.
Survival time = (log <sub>10</sub> labeled population - 2) x labeled D-value	
Kill time / Tiempo de muerte _____	min.
Kill time = (log <sub>10</sub> labeled population + 4) x labeled D-value	
D - value / Valor D (132 °C) _____	sec./seg.
Survival time / Tiempo de sobrevida _____	min.
Survival time = (log <sub>10</sub> labeled population - 2) x labeled D-value	
Kill time / Tiempo de muerte _____	min.
Kill time = (log <sub>10</sub> labeled population + 4) x labeled D-value	
D - value / Valor D (135 °C) _____	sec./seg.
Survival time / Tiempo de sobrevida _____	sec./seg.
Survival time = (log <sub>10</sub> labeled population - 2) x labeled D-value	
Kill time / Tiempo de muerte _____	min.
Kill time = (log <sub>10</sub> labeled population + 4) x labeled D-value	
Z-value / Valor Z _____	°C

Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-3:2017 and IRAM 37102:1999 (Parts 1 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant  
*Terragene®* and *Bionova®* are registered trademarks of Terragene S.A.  
*ATCC®* is a registered trademark of American Type Culture Collection.

<b>Uso exclusivo para profesionales e Instituciones Sanitarias.</b>	
<b>Explanation of Symbols</b>	<span><span></span></span> Manufacture Date
<span><span></span></span> Product is designed for use with Steam sterilization cycles	<span><span></span></span> Expiration Date
<span><span></span></span> Batch number	<span><span></span></span> Manufacturer

### EN Biological Indicators For Steam sterilization

#### Composition

Each tube contains a population of *Geobacillus stearothermophilus* ATCC® 7953 spores inoculated on a carrier. The tube also has a purple growth indicator medium contained within the glass ampoule.

#### Product description

Bionova® BT20 Biological Indicators have been designed for monitoring of vacuum assisted and gravity air displacement Steam sterilization cycles at 121-135 °C. If the sterilization process was not effective, the indicator medium will turn from purple to yellow after incubation between 55-62 °C, thus indicating the presence of living *Geobacillus stearothermophilus* spores. If the sterilization process was effective, the indicator medium will remain purple after incubation. The final readout should be performed at 24 hours of incubation between 55-62 °C.

#### Precautions

Do not use Bionova® BT20 Biological Indicators for monitoring Ethylene Oxide, Dry Heat, Chemical Vapor, Radiation or other sterilization processes.

Do not re-use biological indicators.

Do not use the sterilizer until the biological indicator test results are negative.

#### Instructions for use

1. Identify the Bionova® BT20 Indicator by writing the sterilizer number (in case of having more than one), load number and processing date on the indicator label.

2. Pack the biological indicator along with materials to be sterilized

in an appropriate package according to recommended sterilization practices. Place this package in those areas which a priori you consider most inaccessible for the sterilizing agent (Steam). Generally, the center of the load and areas near the door are problematic.

3. Sterilize as usual.

4. After the sterilization process has finished, open the sterilizer door, wait for 5 minutes and remove the biological indicator from the package. **CAUTION:** Wear safety glasses and gloves when removing the Bionova® BT20 Biological Indicator from the sterilized package. **WARNING:** Do not crush or handle the biological indicator excessively, since this might cause the glass ampoule to burst.

5. Allow the biological indicator to cool down to room temperature.
6. Check the chemical indicator on the label of the biological indicator. A color change to brown confirms that the biological indicator has been exposed to steam. **IMPORTANT:** this color change does not indicate that the process was sufficient to achieve sterility. If the chemical indicator is unchanged, check the sterilization process.

7. Crush the ampoule contained in the biological indicator with an individual ampoule crusher or with the ampoule crusher placed within the top of the incubator s incubation area.

Then shake the tube down vigorously, with movements similar to those performed to lower the temperature of a mercury thermometer, until the medium reaches the base of the tube and soaks the spore carrier entirely. Finally, place the biological indicator in the incubator. **IMPORTANT:** Use a non-sterilized biological indicator as a positive control at least once per day, when a sterilization cycle is run. The positive control ensures that appropriate incubation conditions were met. Both the positive control indicator and the processed indicator should belong to the same batch.

8. Incubate the processed biological indicator and the indicator used as a positive control for a maximum of 24 hours between 55-62 °C. Readout should be carried out at convenient intervals of 10 hours. Color change of the growth indicator medium from purple to yellow means that a sterilization process failure has occurred. If at 24 hours there is no color change in the processed indicators, a final negative result is obtained (the sterilization process was effective). The positive control indicator should show a purple to yellow color change for the results to be valid. **NOTE:** To extend the incubation time beyond 24 hours, we recommend using a humidified environment to prevent complete evaporation of the culture medium contained within the indicator. Record the positive results and discard the biological indicators immediately as described below.

**Storage**  
 Store in a dark place under the next conditions: Temperatures between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

#### Shelf life

Biological indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date. Chemical Process Indicator on SCBI label have an expiration date of 2 years when used as part of SCBI.

Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

#### Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicator can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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