Instructions for Use








---







FORMFIT WALKER & FORMFIT  
WALKER AIR

High Top & Low Top

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	11
IT	Istruzioni per l'uso	14
NO	Bruksanvisning	16
DA	Brugsanvisning	18
SV	Bruksanvisning	20
EL	Οδηγίες Χρήσης	22
FI	Käyttöohjeet	25
NL	Gebruiksaanwijzing	27
PT	Instruções de Utilização	29
PL	Instrukcja użytkownika	31
CS	Návod k použití	34
TR	Kullanım Talimatları	36
RU	Инструкция по использованию	38
JA	取扱説明書	40
ZH	中文说明书	42
KO	사용 설명서	44
RO	Indicații de utilizare	46
SK	Indikácie použitia	48
ET	Kasutusjuhend	50
AR	تعليمات الاستخدام	52

# ENGLISH

---

## SYMBOLS

**MD** Medical Device

## INTENDED USE

The device is intended for immobilization of the foot and ankle

### *Indications for use*

#### **High Top:**

- Foot and ankle conditions that may benefit from motion restriction, such as:
  - Grade 2 and 3 ankle sprains
  - Stable fractures
  - Soft tissue injuries
  - Post-operative, trauma or rehab use

#### **Low Top:**

- Foot and ankle conditions that may benefit from motion restriction, such as:
  - Stable foot fractures
  - Soft tissue injuries
  - Post-operative, trauma or rehab use

No known contraindications.

### **Warnings and Cautions:**

Do not remove the device unless instructed by a healthcare professional.

- Regular healthcare professional supervision is recommended for patients with peripheral vascular disease, neuropathy, and sensitive skin.
- This device should not be applied on patients who are not able to communicate discomfort.
- Patient should use caution when walking on slippery or wet surfaces to avoid injury.
- Ensure the patient performs circulation checks. If a loss of circulation is felt, or if the walker feels too tight, release contact closure straps and adjust to a comfortable level. If discomfort continues, the patient should contact their healthcare professional immediately.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.
- Care should be taken not to overtighten the device.
- Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.

#### **Air version only:**

- Ensure not to overinflate the air chamber.
- Over inflation of air chamber may cause skin compromise in patients with diabetic neuropathy. Deflate air chamber to reduce compression with patients reporting discomfort or showing skin discoloration.
- Ensure patient performs circulation checks. If a loss of circulation is felt, or if the device feels too tight, deflate the air chamber to reduce compression to a comfortable level. If discomfort continues, the patient should contact their healthcare professional immediately.

- In high altitude conditions, the air chamber will expand past its optimal level. Deflate air chamber to reduce compression to a comfortable level.

## GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

## FITTING INSTRUCTIONS

1. Unfasten straps and remove liner from walker (**Fig. 1**). Do not remove pull tabs yet.
2. Place lower leg and foot into the liner and close following the order shown (**Fig. 2**).
3. Insert foot into walker (**Fig. 3**). Ensure that the heel is fully seated in the shell and that the struts are centered over the ankle bone. The foot should be 90° to the lower leg. There should not be a gap between the heel and shell (**Fig. 4**).
4. Place additional pads as needed to accommodate for rubbing and spacing from the liner, shell, or straps (**Fig. 5**).
5. Pull release tabs on both sides on the shell to attach to the liner (**Fig. 6**) and fasten straps as shown (**Fig. 7**).
6. Air versions only: turn the valve clockwise and press the air bladder until the desired inflation is achieved (**Fig. 8, 9**).
7. Remove the walker by unfastening the straps and opening the liner. Air versions only: turn the valve counterclockwise to deflate the liner before doffing.

## Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

## USAGE

### *Cleaning and care*

Device plastic can be washed with a damp cloth and mild soap. Remove liner from the device before washing.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

## DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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