

Instructions for Use

REBOUND[®] AIR WALKER

High top & Low top

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	8
ES	Instrucciones para el uso	10
IT	Istruzioni per l'uso	12
NO	Bruksanvisning	14
DA	Brugsanvisning	16
SV	Bruksanvisning	18
FI	Käyttöohjeet	20
NL	Gebruiksaanwijzing	22
PT	Instruções de Utilização	24
JA	取扱説明書	26
ZH	中文说明书	28
KO	사용 설명서	30

ENGLISH

USES

Rebound Air Walker:

- Foot and ankle conditions that may benefit from motion restriction, such as:
 - Grade 2 and 3 ankle sprains
 - Stable foot fractures
 - Tendon or ligamentous sprain
- During postoperative or rehab phase.
- If used with the Achilles Wedge Kit, Achilles ruptures can be treated as well.

Rebound Air Walker Low Top:

- Foot and ankle conditions that may benefit from motion restriction, such as:
 - Stable foot fractures
 - Tendon or ligamentous sprain
- During postoperative or rehab phase.

CAUTIONS/WARNINGS

READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE USING THE DEVICE.

- **This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.**
- Federal law requires this device to be sold by a physician or by the order of a physician.
- Ensure not to over inflate air chamber.
- Caution should be used if prescribed for patients with diabetes, vascular deficiency, and neuropathy:
 - The skin should be checked regularly
 - Over inflation of air chamber may cause skin compromise in patients with diabetic neuropathy. Deflate air chamber to reduce compression with patients reporting discomfort or showing skin discoloration.
- Ensure patient performs circulation checks. If a loss of circulation is felt, or if the walker feels too tight, deflate air chamber to reduce compression to a comfortable level. If discomfort continues, the patient should contact their medical professional immediately.
- This device should not be applied on patients who are not able to communicate discomfort.
- In high altitude conditions, the air chamber will expand past their optimal level. Deflate air chamber to reduce compression to a comfortable level.
- User should use caution when walking on slippery or wet surfaces to avoid injury.
- User maximum weight limit 300 lbs (200 lbs for the XS version).

PREPARATION

1. Open liner by first unfastening straps and removing front shell.
 - **Note:** Do not remove air bulb from liner.

APPLICATION

2. While seated, place leg in brace.
 - **Note:** Toe Guard can be trimmed off if so desired. (**Fig. 1 & 2**).

- **Note:** Walker circumference (around calf) can be widened by cutting Flex Edge® (**Fig. 3 & 4**)
- 3. Wrap liner securely around foot then leg (**Fig. 5**).
- 4. Place front shell over ankle making sure that the air bulb is aligned with the opening (**Fig. 6**).
- 5. Fasten ankle (Only applies to the high top version), foot, and calf strap in this order. Tighten straps until snug and comfortable.

INFLATION/DEFLATION

6. Inflate air chambers by pressing “+” on the air bulb, and deflate by pressing “-“. Inflate or deflate until snug. (**Fig. 7**)
 - **Caution:** Do Not Over Inflate!
7. Before removing brace, deflate air chambers by pressing “-“.

WASHING INSTRUCTIONS

- Remove liner from walker before washing
- The liner should be washed by hand using a mild detergent and 30°C (86°F) cold water. Do not use bleach or strong detergent. Do not wash in washing machine or place in dryer.
- Line dry in shade.
- Walker plastic can be washed with a damp cloth and mild soap.

LIST OF MATERIALS:

Liner: Polyester, Nylon, Polyurethane, PVC (no bpa), ABS

Shell: Polypropylene, EVA, TPO, Steel, Nylon, Cyanoacrylate

LIMITED WARRANTY

Össur warrants, to the original purchaser, that this product is free from defects in materials and workmanship for six (6) months from initial purchase. This warranty does not apply if this product (1) was not purchased from Össur or its authorized dealer, (2) has been altered in any way, or (3) has not been used in accordance with Össur's Instructions For Use (IFU). This Limited Warranty does not cover damage due to accidents, neglect, misuse, or operation beyond capacity, parts damaged by improper installation, substitution of parts not approved by Össur, or any other alteration or repair by others that, in Össur's judgment, materially or adversely affects the product or part. For complete warranty guidelines please visit www.ossur.com

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