

Diafactory


Dermatophyte Test Strip



Information & Instruction Guide

Diafactory

Dermatophyte Test Strip



- Diagnosis in just 5 minutes
- 97% accurate
- Results unaffected by topical or oral anti-fungal treatments
- Detects 99% of all common dermatophytes with only a 1mg specimen
- No need to send nail sample off to the lab and await results
- Diagnose mycotic infections with confidence and without delay
- Can be used to determine mycological cure

ARTG#: 402694

Will test positive for all common dermatophytes including:

- T. mentagrophytes
- T. rubrum
- T. tonsurans
- T. violaceum
- T. verrucosum
- M. gypseum
- M. canis
- E. floccosum




VIDEO - How to use Diafactory

We have produced the following video to display the simple and quick procedure for the Diafactory® Test Strip:





Diafactory Instructions of Use - Summary

Click here to view full instructions supplied with the Diafactory® Testing Kit:




Each Diafactory® Testing Kit includes the following:

- 1 x Bottle of Extraction Buffer
- 10 x Test Strips
- 10 x Test Tubes
- 10 x Stirring Rods



Diafactory

Diafactory® Test Strips are an immunochromatography kit intended to detect Trichophyton-derived antigens in nails using anti-Trichophyton mouse monoclonal antibody that has been immobilized on a nitrocellulose membrane.



Key Materials

The test strip used in this kit is composed of a sample pad, a reagent strip, a judgment strip and an absorbent pad. The reagent strip contains gold colloid-labeled anti-Trichophyton mouse monoclonal antibody in the dry state, and the judgment strip contains anti-Trichophyton mouse monoclonal antibody in the dry state affixed on the test line zone and the dye in the dry state affixed on the control line zone.

This dye is a colorless dye at a pH of 3 that turns pink at a pH of approximately 4 or higher, and allows the user to confirm that a sample has correctly passed through the test line zone.

From sample to result

A sample that has infiltrated the sample pad (hereinafter called, "the extracted sample") moves to the reagent strip, on which a Trichophyton-derived antigen in the extracted sample binds to a gold colloid-labeled anti-Trichophyton mouse monoclonal antibody to form an immune complex.

While proceeding through the judgment strip, the immune complex is captured by the anti-Trichophyton mouse monoclonal antibody affixed on the test line zone, resulting in the appearance of a purple line of gold colloid (in case it is positive).

If the sample does not contain Trichophyton-derived antigen, no immune complexes are formed and the sample containing unbound gold colloid-labeled anti-Trichophyton mouse monoclonal antibody passes over the test line zone without producing a visible band on the test line zone.

The extracted sample containing unused gold colloid labeled anti-Trichophyton mouse monoclonal antibodies, whether it is Trichophyton-derived antigen positive or negative, passes through the test line zone and reaches the control line zone, where the extracted sample reacts with immobilised dye, resulting in the appearance of a pink band.

Diafactory Research

Click the links below to access research papers on Diafactory®:

ORIGINAL ARTICLE
Clinical study of Dermatophyte Test Strip, an immunochromatographic method, to detect tinea unguium dermatophytes
 Yuchiko TSUNEMI¹, Masatosi HIRUMA²
¹Department of Dermatology, Tokyo Women's Medical University, ²Chunichi Institute for Medical Mycology and Allergology, Tokyo, Japan

ABSTRACT
 The Dermatophyte Test Strip visualizes mycotic antigens by immunochromatography. It allows easy and fast detection of dermatophytes. A multicenter, single-arm, comparative clinical study was designed to evaluate the capacity of Dermatophyte Test Strip to detect dermatophytes in suspected tinea unguium specimens in comparison with direct microscopy and polymerase chain reaction (PCR). Signed consent was obtained from 222 subjects and all subjects completed the study. With the Dermatophyte Test Strip, dermatophytes were detected in 201 of 222 (90.5%) specimens but not in 22 of 222 (9.5%) specimens. With direct microscopy, dermatophytes were detected in 170 of 222 (76.6%) specimens but not in 52 of 222 (23.4%). Of the 40 specimens that showed inconclusive results between the two methods, PCR gave further results for 40 specimens, of which 27 (67.5%) specimens were positive and three (7.5%) were negative for dermatophytes. The positive concordance rate, negative concordance rate and overall concordance rate between the Dermatophyte Test Strip and direct microscopy were 81.1%, 68.7% and 79.1%, respectively. When inconclusive results were corrected using the results of PCR, these rates were 82.0%, 71.4% and 86.0%, respectively. When the specimens that could not be tested by PCR because no pieces for the PCR test was left were excluded from analysis, these rates were 99.2%, 78.0% and 97.2%, respectively. The present results indicate good detection capacity of the Dermatophyte Test Strip. The Dermatophyte Test Strip provides a reliable, convenient and quick method to test for tinea unguium.

Key words: dermatophytes, immunochromatography, immunological diagnosis, tinea unguium, Trichophyton.

INTRODUCTION
 Tinea unguium is a nail disease caused by dermatophyte infection of the nail plate.¹⁻³ Its prevalence in advanced countries is considered as at least 10% of the population,⁴ and is about 1% in areas of approximately 10%.⁵ The incidence increases with aging in all countries, and it is therefore considered to be a disease of aging. Patients with tinea unguium pose a problem of recognizing the source of infection of their own feet. Tinea unguium is particularly considered as a predictor of diabetic foot syndrome.⁶ Thus, tinea unguium does not only cause problems with nail appearance, but it also seriously impairs patients' quality of life. It should therefore be treated actively as far as possible.⁷

Because many other diseases have similar symptoms as tinea unguium,⁸ differentiation is not easy, and a definite diagnosis by mycological examination is necessary before treatment is started. Dermatophytes are classified into three genera, namely, *Trichophyton* (T), *Microsporum* (M) and *Epididymomyces* (E).⁹ T culture accounts for the primary cause of tinea unguium,¹⁰ followed by T and M.¹¹

The diagnostic tests for tinea unguium include direct microscopy with potassium hydroxide (KOH), fungal culture, periodic acid-Schiff staining (PAS), fluorescence staining, confocal microscopy and molecular biological methods such as polymerase chain reaction (PCR).¹² Although direct microscopy with KOH and fungal culture have been the gold standard for the diagnosis of tinea unguium,¹³ these tests sometimes are a burden in the clinical setting, because direct microscopy requires experience to identify the fungal elements, and diagnosis of tinea in time-consuming. It takes 2-3 weeks to obtain results using fungal culture, and the detection rate is lower than with direct microscopy.¹⁴ Therefore, the development of an easy and quick test to diagnose tinea unguium accurately is a high demand.

The recently developed Dermatophyte Test Strip detects dermatophytes easily and rapidly by immunochromatography.

Correspondence: Yuchiko Tsunemi, M.D., Ph.D., Department of Dermatology, Tokyo Women's Medical University, 8-1 Kawasumi, Shinjuku-ku, Tokyo 162-8601, Japan. Email: ytsunemi@twmu.ac.jp
 Received 13 October 2015; accepted 15 January 2016
 © 2016 Japanese Dermatological Association 1417

Clinical study of Dermatophyte Test Strip, an immunochromatographic method, to detect tinea unguium dermatophytes

Journal of Medical Diagnostic Methods
 Volume 5 Issue 2 • 1001216
 DOI: 10.4137/JMD.M1001216

Research Article
Development of a New Dermatophyte-Detection Device using Immunochromatography
 Hiroaki Watanabe¹ and Masatosi Miyama²
¹Life Chemistry Research Office, JNC Corporation, Tokyo, Japan
²Chunichi Institute for Medical Mycology and Allergology, Life Chemistry Research Office, JNC Corporation, Shin-Okubo Building, 2-1, Chitose, Chitose, Chiyoda-ku, Tokyo 1001021, Japan. Tel: 03-6342-0000. E-mail: h.watanabe@jnc.com

Received 26 May 2016; Accepted 04 Jun 2016; Published 04 Jun 2016
 Copyright © 2016 Watanabe et al. This is an open access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Abstract
 Background: The traditional method for the diagnosis of tinea unguium is direct microscopy. However, with anti-dermatophyte monoclonal antibody have been developed that recognize the cell wall polysaccharide antigen of dermatophyte fungi including dermatophytes, we applied it to immunochromatography. We describe the development of the technology, proposed it, and demonstrate that this test method provided useful information for a diagnosis of tinea unguium by analyzing large numbers of clinical specimens. However, these studies were carried out by a dermatologist.

Objective: Establishment of the usage of the test strip for clinical use and the certification of the stability during long-term storage.

Methods: Various fungi and bacteria were cultured and extracted by extraction buffer to evaluate the reactivity of monoclonal antibody of dermatophytes. Trichophyton solitum was used as a quality control antigen, and detection limits were set at 0.1 µg/ml as used positive and set at a 100-fold concentration (100 µg/ml) as a negative control. Test antigens, which had already been identified as positive or negative using the test strips, were randomly selected, cut into five pieces and mixed. Both positive and negative standard test samples were prepared.

Results: Positive or negative results were obtained from four test strips at 30°C and 4°C. The detection limits of dried 7 samples of dermatophytes were 0.1 to 1 µg/ml. The results of reactivity showed that the detection limits were positive. On the other hand, the test strips did not react to *Microsporum* or *Candida* species and bacterial strains. Some of *Aspergillus*, *Penicillium*, and *Fusarium*, which are usually not resistant microorganisms that grow in healthy humans, showed a positive reaction. The antifungal agents, miconazole, griseofulvin, and itraconazole did not affect the results. All lots of the test strips met the standards by the method of quality control after they were kept in storage at 30°C for up to 22 months.

Conclusion: A newly developed dermatophyte-detection device was easy to use, gave rapid results and high reactivity, and was stable for 22 months at 30°C.

Keywords: Dermatophytes; Monoclonal antibody; Device

Introduction
 Tinea unguium is a nail disease caused by dermatophytes [1,2]. The traditional method for a diagnosis is mainly KOH direct microscopy and fungal culture, which comprise the gold standard. It takes 2 to 3 weeks for fungal culture and the results are obtained, and the accuracy of the culture is approximately 20-70%, which is considered to be although low [3]. Therefore, microscopy is usually used to diagnose tinea unguium in clinical practice. However, a rigorous microscope and the operator's experience is necessary for microscopy, and there is also the problem that it takes a long time to identify the test [4,5].

Nishi et al. have developed an anti-dermatophyte monoclonal antibody that recognizes the cell wall polysaccharide antigen of dermatophyte fungi, including dermatophyte, and applied it to immunochromatography. They then subjected 37 nail specimens to the test and successfully detected dermatophyte antigen [6,7]. Tamami et al. also reported that the test strip could be used to identify the test with large numbers of clinical specimens, and also demonstrated that the test method provided useful information for a diagnosis of tinea unguium [8].

However, these studies were carried out by a sensitive procedure of the test strip, and to put this test strip into clinical practice widely as an in vitro diagnostic, it was necessary to determine a practical procedure of use along with the issue of cost and clearly understand the reproducibility of the results, and long-term storage stability of the test strip.

In this report, we established the usage of this test strip for practical application, and investigated the monoclonal antibody and long-term storage stability.

Volume 5 Issue 2 • 1001216
 11 May 2016
 ISSN 2258-9754 (print), an open access journal

Development of a New Dermatophyte Detection Device using Immunochromatography

CLINICAL AND LABORATORY INVESTIGATIONS
 British Journal of Dermatology

Screening for tinea unguium by Dermatophyte Test Strip*
 Y. Tsunemi¹, Y. Takekura², Y. Miura², K. Nakagami³, H. Sanada³ and M. Kawakami³
¹Department of Dermatology, Tokyo Women's Medical University, 8-1 Kawasumi, Shinjuku-ku, Tokyo 162-8601, Japan
²Department of Ophthalmology, Niigata City University, Niigata City, Niigata, Japan
³Department of Dermatology, Niigata City University, Niigata City, Niigata, Japan

Summary
 Background: The direct microscopy, fungal culture and histopathology that are necessary for the definitive diagnosis of tinea unguium are disadvantages in that detection sensitivity is affected by the level of skill of the person who performs the testing, and the procedures take a long time.
 Objective: The Dermatophyte Test Strip, which was developed recently, can easily and easily detect *Trichophyton* fungi in samples in a short time, and does not require special care for its use as a method for tinea unguium screening. With this in mind, we examined the detection capacity of the Dermatophyte Test Strip for tinea unguium. Methods: The presence or absence of fungal elements was judged by direct microscopy and Dermatophyte Test Strip in 145 nail samples obtained from outpatients in screening rooms for the elderly. Moreover, the minimum sample amount required for positive determination was estimated using 33 samples that showed positive results by Dermatophyte Test Strip.
 Results: The Dermatophyte Test Strip showed an 89% sensitivity, 10% specificity, 84.8% positive predictive value, 97% negative predictive value and a positive and negative concordance rate of 89.1%. The minimum sample amount required for positive determination was 0.0025–0.212 mg.
 Conclusion: The Dermatophyte Test Strip showed very high sensitivity and negative predictive value, and was considered a potentially useful method for tinea unguium screening. Positive determination was considered to be possible with a sample amount of about 1 mg.

What's already known about this topic?
 • Direct microscopy, culture and histopathology are employed to detect fungal diseases.
 • These methods require technical skill and/or consume time.

What does this study add?
 • The Dermatophyte Test Strip visualizes mycotic antigens by immunochromatography.
 • It allows easier and faster detection of fungi in samples, with very high sensitivity.
 • It is a useful method for screening of tinea unguium.

Tinea unguium is estimated to occur in at least 10% of the population in advanced countries,¹ and is increasing along with aging in all countries. Thus, tinea unguium is considered to be a disease of aging. Patients with tinea unguium pose a problem of recognizing the source of infection of their own feet. Tinea unguium is particularly considered as a predictor of diabetic foot syndrome.² Thus, tinea unguium does not only cause problems with nail appearance, but it also seriously impairs patients' quality of life. It should therefore be treated actively as far as possible.³

Because many other diseases have similar symptoms as tinea unguium,⁴ differentiation is not easy, and a definite diagnosis by mycological examination (direct microscopy, fungal culture and histopathology) is necessary in order to start the treatment of tinea unguium. However, there are disadvantages in direct microscopy, such as detection sensitivity being affected by the level of skill of the person who performs the test. Fungal culture takes 2–3 weeks

DOI: 10.1111/bjd.12660
 © 2016 British Association of Dermatologists

Screening for tinea unguium by Dermatophyte Test Strip

GUIDELINES
 British Journal of Dermatology

British Association of Dermatologists' guidelines for the management of onychomycosis 2014
 M. Bennett,^{1,2} J. Lear,³ V. Mahida,⁴ R. H. M. Hoogkamp⁵ and M. Richardson^{6,7}
¹Nuffield Health NHS Foundation Trust, Nail Clinic, London N9F 3DQ, U.K.
²Manchester Academic Health Science Centre (MAHSC), University of Manchester, 46 Castle Street, Manchester M1 1WT, U.K.
³Galderma, Galderma (UK) Healthcare, Ltd, London, London W6 0LE, U.K.
⁴British Association of Dermatologists, Wilton House, 4 Foyers Square, London W1T 5HQ, U.K.
⁵Myking Medical Care, Manchester, University Hospital of South Manchester, Salford Road, Manchester M13 9PL, U.K.

1.0 Purpose and scope
 The overall objective of the guidelines is to provide up-to-date, evidence-based recommendations for the management of onychomycosis. The document aims to (i) offer an overview of all relevant literature since January 2003; focusing on any key developments; (ii) address important practical clinical questions relating to the primary guideline objectives, for example accurate diagnosis and identification of cases and suitable treatment to minimize the duration of disease and (iii) provide guideline recommendations and, where appropriate, with some health economic implications; and (iv) discuss potential developments and future directions.

The guidelines are presented as a detailed review with high-level recommendations for practical use in the clinic, in addition to an updated patient information leaflet (available on the British Association of Dermatologists' (BAD) website, www.bad.org.uk).

2.0 Stakeholder involvement and peer review
 The guideline development group consisted of consultant dermatologists and a consultant mycologist. The draft document was circulated to the BAD membership, the British Dermatological Nursing Group, the Primary Care Dermatological Society and the North West Region Kirby Foster Association for comment, and was peer reviewed by the Clinical Standards Unit of the BAD (lead up of the Therapy & Guidelines Subcommittee) prior to publication.

3.0 Methodology
 The set of guidelines has been developed using the BAD's recommended methodology⁸ and with reference to the Approval of Guidelines: Research and Evaluation (AGREE II) instrument (www.agreement.org). Recommendations were developed for implementation in the National Health Service using a process of considered judgement based on the evidence. The Preferred, Modified and Further discussion were searched for (see online, additional and non-recommended considered clinical trials, case reports, case series and open studies involving onychomycosis published in the English language from January 2003 to the present).

DOI: 10.1111/bjd.12158
 © 2014 British Association of Dermatologists

British Association of Dermatologists' guidelines for the management of onychomycosis 2014

Diafactory

Dermatophyte Test Strip

briggatemedicalcompany

Ph: (03) 8586 7800
Free Call: 1800 33 4142
Fax: (03) 8586 7888

briggate@briggate.com.au

www.briggatemedical.com