



Life Span:
80,000 Hours

Warranty:
5 years



AFMA L.E.D Magnifying Lamp & Stand

AFMAPACK-2


Afma LED Bulk Package

1 x Afma LED Lamp (in White or Silver Grey)

1 x Afma Telescopic Mobile Stand

- Brighter, whiter light resulting in clearer subject definition
- Energy efficient - uses less power than fluoro tubes
- Low on-going costs - no tube or starter switch replacement
- Long lasting - *Average life span of 80,000 hours*
- SAFE:
 - operates at a low temperature
 - mercury, UV & IR free
- 3 diopter lens
- Handle design light housing for greater positioning control
- 5 Year Warranty on AFMA L.E.D Lamp

Log-in to www.briggatemedical.com for pricing



| Mobile Stand Dimensions: | |
|--------------------------|--------|
| Base Width | 600mm |
| Max. Height | 1350mm |
| Min. Height | 880mm |



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com