



Manufactured in  
GERMANY

A perfect product is not one person's idea, but rather the result of constant respectful collaboration. Professional qualifications and technical know-how form the foundation. Passion for innovative research and personal motivation to be among the best – it is only this that makes us a successful global player in precision rotational instruments. The highly sensitive industry in which we move demands our maximum care at every stage of development. Our employees carry this almost painstaking basic attitude from research through extensive test phases through to the production of each individual product with conviction. They are very aware of their responsibility towards our customers and their patients.

Progressive reliability and cautious growth in a motivating, collegial atmosphere – characteristics of Acurata, which form the foundation for our more than 35-year record of success.

Karl Mahnhardt  
Managing Director  
Acurata GmbH

### ***Not all instruments are created equal...***

Podiatry drills require tools with a shaft thickness of approx. 2.35mm. It is also important that the tool shafts are as straight as possible. Over time, tool thickness and straightness can be effected, and thus the chances of damaging the clutch mechanism of you drill are increased. Regular testing of your tools to meet these requirements can be achieved with the use of the Acurata Concentricity Tester.

