

# AC Blue Tungsten

AC Blue is a specifically formulated coating designed to protect and enhance the burr. Using a process called Physical Vapour Deposition (PVD), a micron-thin layer of AC Blue is permanently bonded to each workpiece, resulting in:



**Greater durability – longer service life**



**Reduced heat build-up**



**Easier cleaning**

- Allows efficient debridement of mycotic and thickened nails
- Creates less dust particles
- Long lasting burr made from hardened tungsten carbide alloy
- Cutting edge remains sharp



## MEDIUM CUT



**BTBUDM60B**  
Medium Cut  
Medium Bud



**BTCONEM60B**  
Medium Cut  
Cone

## COARSE CUT



**BTBUDC60B**  
Coarse Cut  
Medium Bud



**BTCONEC60B**  
Coarse Cut  
Cone

# briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987