

Ceramic Burrs

- Designed to leave a smooth, clean surface after completion of treatment
- Easy to clean and autoclavable
- Offers reduced heat build-up and reduced vibration

CROSS CUT

Cross Cut offers high removal rates of callus and shaping of mycotic nails



BC1907460
Cross Cut
Medium Cone




BC1909440
Cross Cut
Narrow Cone




COARSE CUT

Coarse Cut offers VERY high removal rates of callus and thickened nails



BCLSQ7460
Coarse Cut
Medium Cone



BCLSQ9440
Coarse Cut
Narrow Cone



BCLSQ1060
Coarse Cut
Cylinder



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987