

# Hurricane Burrs

- Enhances the dust control capabilities of your drill
- Fast and effective callus removal
- Long lasting diamond grit coating
- Alternative to Sanding Discs

## COARSE GRIT



**BDHURS**  
Coarse Grit  
Small



**BDHURM**  
Coarse Grit  
Medium



**BDHURL**  
Coarse Grit  
Large




## SUPER COARSE GRIT



**BDHURSCM**  
Super Coarse Grit  
Medium



**BDHURSL**  
Super Coarse Grit  
Large



Designed and constructed with ventilation holes to create a cyclonic extraction effect that enhances the collection of debris. Also aids in the prevention of heat build-up.

# briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987