

BAEHR®

BAEHR TEC A2000+

briggatedrills
SALES • SERVICE • SUPPORT

PERFORMANCE GUARANTEE

contact us for further information



Baehr Tec A2000+ Podiatry Drill with Dust Extraction

- 40,000 rpm Handpiece
- Handpiece On/Off Switch
- LCD Display
- Clockwise/Counter-clockwise Rotation
- Suction Power up to 230L/min
- Dual Brushless Vacuum System
- Advanced Filtration System
- 3 Programmable Memory Buttons for Suction Level, Burr Speed & Burr Direction
- Integrated Handpiece Holder

- Height 10.5cm
- Width 28cm (includes h/p holder)
- Depth 24cm
- Weight 2.1kg
- Handpiece Weight 100g

BAEH-A2K+

Baehr Tec A2000+ Podiatry Drill

Log-in to www.briggatemedical.com for pricing

Suction Performance:
260 L/min

“Silent Suction”

Advanced 5-layer
Filtration System

Warranty:
2 years

Manufactured in:
Germany



BAEHR®
ADVANCED
FILTRATION
SYSTEM

5-layer cartridge system that actively captures dust and bacteria particles and controls odour compounds.



BAEHR®
BRUSHLESS
NEXT GENERATION
MINI TURBINE

x2

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com