

BAEHR®

BAEHR TEC A700

briggatedrills
SALES • SERVICE • SUPPORT

PERFORMANCE
GUARANTEE

contact us for further information



Baehr Tec A700 Podiatry Drill with Dust Extraction

- 27,000 rpm Handpiece
 - Handpiece On/Off Switch
 - LCD Display
 - Clockwise/Counter-clockwise Rotation
 - Suction Power up to 210L/min
 - Brushless Vacuum Motor
 - Multi-Layer Filtration System
 - 2 Programmable Memory Buttons for Suction Level, Burr Speed & Burr Direction
 - Whisper quiet operation: 56dB
 - Easy-to-clean fool-proof chuck
-
- Height 10cm
 - Width 26cm
 - Depth 18cm
 - Weight 1.8kg
 - Handpiece Weight 102g

BAEH-A700

Baehr Tec A700 Vacuum Drill

Log-in to www.briggatemedical.com for pricing

Suction Performance:
210 L/min

“Silent Suction”

Advanced 4-layer
Filtration System

Warranty:
2 years

Manufactured in:
Germany



BAEHR®
ADVANCED
FILTRATION
SYSTEM

4 layer cartridge system that actively captures dust and bacteria particles and controls odour compounds.



BAEHR®
BRUSHLESS
NEXT GENERATION
MINI TURBINE

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com