



Includes zip-up
"sports style" carry bag



Baehr Mobile Legrest with Carry Bag

- Minimum Height: 510mm
- Maximum Height: 760mm
- Legrest Pad Dimensions: 210mm x 150mm
- Unit Weight: 1.3kg
- Safe Working Load: 12.8kg

- Aluminium frame construction
- Angle adjustable legrest pad
- Rotary lever adjusts leg rest height
- Foldable tripod style
- Includes Carry Bag

LEGRESTBAE

Baehr Mobile Legrest with Carry Bag

Log-in to www.briggatemedical.com for pricing

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987