

# BERCHTOLD. podoQ 40



## Berchtold PodoQ 40 Podiatry Drill with Dust Extraction

- 25,000 rpm Handpiece
  - Handpiece On/Off Switch
  - LCD Touch Screen Display
  - Clockwise/Counter-clockwise Rotation
  - Brushless Vacuum Motor Technology
  - Touch Button Controls
  - 3 Programmable Memory Buttons for Suction Level, Burr Speed & Burr Direction
  - Handpiece Holder
- 
- Height 8cm
  - Width 33cm (includes h/p holder)
  - Depth 25cm
  - Weight 2.4kg
  - Handpiece Weight 140g

### BERCH-Q40

Berchtold PodoQ 40 Podiatry Drill  
with Dust Extraction

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing

Suction Performance:

190 L/min

Warranty:

3 years

Manufactured in:

Germany



**briggatedrills**  
SALES • SERVICE • SUPPORT

**PERFORMANCE  
GUARANTEE**

contact us for further information

# briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

[briggate@briggate.com.au](mailto:briggate@briggate.com.au)

[www.briggatemedical.com](http://www.briggatemedical.com)