

Hadeco Bidop 7

BI-DIRECTIONAL DOPPLER

Warranty:
1 year

Manufactured in:
Japan



LARGE COLOUR SCREEN

89mm colour LCD screen shows high resolution real-time waveforms and numerical data

BIGGER BATTERY

Larger capacity Ni-MH battery for longer operational time

USER-FRIENDLY

Improved menus and button configuration for ease-of-use

IMPROVED CONNECTIONS

Snap-lock probe connectors make swapping probes quick and easy

IMPROVED AUDIO

Larger speaker with improved technology for clean and crisp waveform audio

WIDE BEAM 8MHZ PROBE

Higher sensitivity wide beam probe enhances audio quality

BLOOD PRESSURE MODE

Easier PPG and arterial waveform interpretation



- LARGER 89mm COLOUR LCD Real-time Waveform Display
- HIGHER sensitivity wide beam 8MHz Probe
- IMPROVED speaker technology for high definition audio
- LARGER capacity Ni-MH Battery with Charger
- User friendly menu system and navigational buttons
- 30 Waveform Memory
- Snap-lock probe connector for easy probe swapping
- Blood volume flow calculations

Dimensions:

Width 11.5cm
Depth 21.5cm
Height 6.5cm
Weight 570g


Compatible with:

- PPG Probe
- Smartlink Software
- Temperature Probe

DOP-BIDOP7

Log-in to www.briggatemedical.com for pricing

ADD
SMART-V-LINK
SOFTWARE WHEN
PURCHASING THE
HADECO BIDOP 7 &
SAVE \$75.00+GST



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com