

# Hadeco Bidop 7 T.B.I Kit

DOPPLER



SPHYGMO



Contains everything required to obtain ABI/TBI readings:

- 1 x BiDop 7 Bi-Directional Doppler
- 1 x Hadeco PG30 PPG Probe with snap-lock connectors
- 1 x Riester Ri-san Sphygmomanometer with Adult Cuff (24cm-34cm)
- 1 x 2.5cm Digital Cuff (latex-free)
- 1 x Spring Clip to secure PPG Probe to Toe


#### Hadeco Bidop 7 Specifications:

- LARGER 89mm COLOUR LCD Real-time Waveform Display
- HIGHER sensitivity wide beam 8MHz Probe
- IMPROVED speaker technology for high definition audio
- LARGER capacity Ni-MH Battery with Charger
- User friendly menu system and navigational buttons
- 30 Waveform Memory
- Snap-lock probe connector for easy probe swapping
- Blood volume flow calculations

#### DOP/PACK7

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing

**ADD SMART-V-LINK  
SOFTWARE FOR  
\$275.00+GST  
(NORMALLY \$350.00+GST)**



# briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

[briggate@briggate.com.au](mailto:briggate@briggate.com.au)

[www.briggatemedical.com](http://www.briggatemedical.com)