

# Carbon Fibre Foot Plates - CONTOURED

## Carbon Fibre Foot Plates

- High strength carbon fibre composite foot plates used to limit forefoot motion
- Suitable for conditions such as hallux rigidus, turf toe, plantar plate pathologies, 1st MPJ arthritis or capsulitis
- Thermoplastic composite laminate made using a semi-crystalline polyethylene terephthalate (PET) matrix resin
- Highly resilient - will not lose it's shape or material properties after continued use in high loading applications
- Low profile design for easy shoe fit
- Minor grinding can be utilised to take back device width
- Sold individually - NOT IN PAIRS



## Carbon Fibre Foot Plates - CONTOURED


\* fits UP TO shoe size shown \* approx. foot plate length

### SEMI-RIGID (Thickness: 1.69mm - 1.91mm)

CFC23-SR	♂ *37/*36	US*: 5	♀ 6	23cm*
CFC26-SR	♂ *41/*40	US*: 8	♀ 9	25.6cm*
CFC29-SR	♂ *44½/*45	US*: 11	♀ 12	28.3cm*
CFC31-SR	♂ *47/*48	US*: 13	♀ 14	30.8cm*

### RIGID (Thickness: 1.93mm - 2.19mm)

CFC23-R	♂ *37/*36	US*: 5	♀ 6	23cm*
CFC26-R	♂ *41/*40	US*: 8	♀ 9	25.6cm*
CFC29-R	♂ *44½/*45	US*: 11	♀ 12	28.3cm*
CFC31-R	♂ *47/*48	US*: 13	♀ 14	30.8cm*



Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing

**briggatemedical**company

Supporting Podiatrists with integrity and trust since 1987