

Carbon Fibre Foot Plates - MORTON'S EXT.

Carbon Fibre Foot Plates

- High strength carbon fibre composite foot plates used to limit forefoot motion
- Suitable for conditions such as hallux rigidus, turf toe, plantar plate pathologies, 1st MPJ arthritis or capsulitis
- Thermoplastic composite laminate made using a semi-crystalline polyethylene terephthalate (PET) matrix resin
- Highly resilient - will not lose its shape or material properties after continued use in high loading applications
- Low profile design for easy shoe fit
- Minor grinding can be utilised to take back device width
- Sold individually - NOT IN PAIRS




Carbon Fibre Foot Plates - MORTON'S

* approx. foot plate length

MORTON'S FLAT

CFMF-SR20	SEMI-RIGID	FLAT	20.2cm*
CFMF-SR26	SEMI-RIGID	FLAT	26cm*
CFMF-R20	RIGID	FLAT	20.2cm*
CFMF-R26	RIGID	FLAT	26cm*




FLAT

MORTON'S CONTOURED - LEFT

CFMCL-R20	SEMI-RIGID	CONTOURED	20.2cm*
CFMCL-R26	SEMI-RIGID	CONTOURED	26cm*
CFMCL-SR20	RIGID	CONTOURED	20.2cm*
CFMCL-SR26	RIGID	CONTOURED	26cm*

MORTON'S CONTOURED - RIGHT

CFMCR-SR20	SEMI-RIGID	CONTOURED	20.2cm*
CFMCR-SR26	SEMI-RIGID	CONTOURED	26cm*
CFMCR-R20	RIGID	CONTOURED	20.2cm*
CFMCR-R26	RIGID	CONTOURED	26cm*



CONTOURED

Log-in to www.briggatemedical.com for pricing

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com