

## WHY IT WILL HELP?

FOR LATERAL FOOT, ANKLE AND PERONEAL PATHOLOGIES

### Stable Rearfoot Control

Provides stability for transition from heel strike to midstance during foot loading

### Lateral Arch Cushioning and Contour

Provides improved shock absorption, stabilises the foot and promotes medial transition to propulsion

### Lateral Wedge

Improves lateral stability and medial transition to propulsion

### Medial Arch Support

Supports the transition from midstance into propulsion, improves windlass activation, forefoot stability and high gear toe off

### Moderate Heel Raise

Reduce tissue stress on posterior structures and improves sagittal plane transition from midstance to propulsion

### Full Length Performance EVA Orthotic

Provides optimal cushioning and control

### Durafit Top Cover

Luxurious, breathable and durable

