

PROPULSION

WHY IT WILL HELP?

IMPROVES ANKLE AND FOREFOOT SAGITTAL PLANE FUNCTION

Stable Rearfoot Control

Provides stability for transition from heel strike to midstance during foot loading

Medial Arch Support

Supports the transition from midstance into propulsion, improves windlass activation, forefoot stability and high gear toe off

Lateral Arch Support

Stabilises the foot and aids the transition from heel strike to toe off

Moderate Heel Raise

Reduce tissue stress on posterior structures and improves sagittal plane transition from midstance to propulsion

Reverse Mortons Extension

Improves propulsion, activation of windlass mechanism and sagittal plane function of the foot

Full Length Performance EVA Orthotic

Provides optimal cushioning and control

Durafit Top Cover

Luxurious, breathable and durable

