

WCS® Light Closed Wound Care Shoe System



WOUND CARE
CIRCLE

- > for treatment of the diabetic foot syndrome
- > for healing of open wounds and minor ulcerations
- > for pressure relief
- > wide toe box for any remaining bandages



DARCO

Pressure points and foot ulcers occur most frequently with diabetes mellitus and in many ways present a great challenge in healing and off-loading sensitive feet.



- > seamless toe area
- > elastic and breathable upper
- > extra space for light bandages in toe area
- > 4 removable insoles in different colours and shore grades; the individual customization of the insoles enables a selective creation of plantar off-loading zones
- > comfortable fitting provided by high heel cup
- > decent look




WCS® Closed Wound Care Shoe System

The diabetic patient requires specifically at an early stage of diabetic foot syndrome with minor ulceration a comfortable and reliable shoe ensuring effective pressure relief for indoor and outdoor use.

The new WCS® Light of DARCO offers the proven four-layer-insole-system with individual customization for off-loading included in a closed shoe. The ultra soft and breathable upper provides maximum comfort and leaves extra space for bandages.

Sizes and Measurements

All measurements are defined in cm and correspond to the insole. DARCO shoe sizes represent a general reference only. The choice of shoe size is significantly influenced by the volume of the bandages and size of cast.



The WCS® is available pairwise in seven sizes.

ART. NO.	COLOR	SIZE	
WCSSL0	Black	XS	35.0 – 36.0
WCSSL1	Black	S	37.0 – 38.0
WCSSL2	Black	S/M	39.0 – 40.0
WCSSL3	Black	M	41.0 – 42.0
WCSSL4	Black	L	43.0 – 44.0
WCSSL5	Black	XL	45.0 – 46.0
		XXL	47.0 – 48.0