

Hadeco ES100VX

Warranty:
1 year

Manufactured in:
Japan



- Unidirectional Hand-Held Doppler
- 8MHz Probe
- Convenient Probe Button freezes data for notation
- Speaker Output: 300 mW +
- Supplied with
 - Headphones
 - Carry Bag
 - 1 x 9V Battery

DOP-MINI

Log-in to www.briggatemedical.com for pricing

Dimensions:

Width	7.8cm
Depth	14cm
Height	2.7cm
Weight	240g



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com