

# THE EQUINUS BRACE

Effective solution for Equinus related pathologies



Significantly Improve Patient Compliance

Engage windlass with 60° Toe Wedge to lock Rearfoot

Lock the Knee into Full Extension to stretch Gastrocnemius

1 Hour a Day, 12 Week Treatment

Adjustable Ankle Dorsiflexion Angle

Available in Adult & Paediatric Sizing

ADULT BRACE      PAEDIATRIC BRACE

## Equinus Related Pathologies and Literature



Plantar Fasciitis



Achilles Tendinopathies



Medial Tibial Stress Syndrome



Metatarsalgia



Sever's Disease




Adult & Paediatric Flatfoot

**EXTENDS PROXIMAL TO KNEE JOINT TO INCORPORATE GASTROCNEMIUS STRETCH.**



**STRETCHES THE WHOLE LOWER LIMB POSTERIOR CHAIN.**



### Equinus Brace

- Simple, effective brace for relief from Equinus related pathologies
- Offers great patient compliance - use 1 hour a day
- Easy-to-use - can be fitted within minutes
- Size adjustable to suit patient limb length

**EQBRACEADU** Adult Size

**EQBRACEPED** Paediatric Size

[CLICK HERE FOR DEMONSTRATIONAL VIDEOS & LITERATURE](#)

**briggate**medicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com