

Hadeco

ES100V3 T.B.I. Kit



Contains everything required to obtain ABI/TBI readings:

- 1 x Hadeco ES100V3 Doppler
- 1 x Hadeco PG21 PPG Probe
- 1 x Riester Ri-san Sphygmomanometer with Adult Cuff (24cm-34cm)
- 1 x Latex-Free Digital Cuff (2.5cm)
- 1 x Spring Clip to secure PPG Probe to Toe




Hadeco ES100V3 Specifications:

- Bidirectional Hand-Held Doppler
- Backlit LCD Screen displays realtime waveforms and numerical data
- 8MHz Probe
- Convenient Probe Button freezes data for notation
- 30 Waveform Memory
- 9V alkaline battery
- Speaker Output: 300 mW +
- Blood volume flow calculations

DOP/PACKV3

Log-in to www.briggatemedical.com for pricing

**ADD SMART-V-LINK
SOFTWARE FOR
\$275.00+GST
(NORMALLY \$350.00+GST)**



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com