


# Flow MINI



## Flow MINI Percussion Therapy Device

MINI gives Athletes a Precision Sports Massage via four levels of Powerful Vibration – to promote best Performance, Reduce Fatigue, and Relieve Muscle Swelling and Tension.

- Accelerates warmup of muscles before exercise.
- Regular use improves an athlete's overall mobility and range of movement.
- Helps to release trigger points and assists the break up of scar tissue.
- Easy-to-use self-myofascial release any time, anywhere.
- May help reduce common muscle and joint pain

### Top Features

- Next Generation ultra quiet high-torque motor (60-65dbs).
- Easy to use 4-speed control with Flow Wave Mode.
- 4 Flow Heads™ included – Medium, Hammer, Bullet & Fork.
- Long battery life and fast recharge time.
- Highly portable and Lightweight design.
- Reliable, robust and hard wearing construction

### Specifications

- Battery life: 2-3 hours with typical usage
- Recharge time: 3 hours
- Charging Method: A/C Wall Charger
- Battery Type: 1500 Mah Rechargeable Lithium-ion
- Level 1 mode: 1800 percussions per minute 'PPM'
- Level 2 mode: 2400 PPM
- Level 3 Mode: 3000 PPM
- Flow mode: 1800 through 3000 PPM

### FLOWMINI

Flow MINI Percussion Therapy Device

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing

A FLOW HEAD FOR EVERY NEED.



#### BALL HEAD

With forging high-density foam, this head is for use around bones and other sensitive areas.



#### HAMMER HEAD

The jack of all trades head, this head is designed for overall body treatment and larger muscles groups.



#### BULLET HEAD

The small diameter of the bullet head makes it perfect for focused muscle treatment; to pinpoint trigger points and work them away.



#### FORK HEAD

Perfect for the spine, lower neck and Achilles' tendon along with major muscle groups like the calves & quadriceps.

# briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

[briggate@briggate.com.au](mailto:briggate@briggate.com.au)

[www.briggatemedical.com](http://www.briggatemedical.com)