



Heel Raisers

Used to reduce strain on the Achilles Tendon (Achilles Tendonitis) or for limb length discrepancies.
Available in heights of 4mm and 6mm (measurement taken from the centre of the heel).
6mm is the most commonly used height. Self adhesive.

FTRAIS4S	Red	Small	4mm Low	Pkt 5 Pairs
FTRAIS4M	Red	Medium	4mm Low	Pkt 5 Pairs
FTRAIS4L	Red	Large	4mm Low	Pkt 5 Pairs
FTRAIS6S	Red	Small	6mm High	Pkt 5 Pairs
FTRAIS6M	Red	Medium	6mm High	Pkt 5 Pairs
FTRAIS6L	Red	Large	6mm High	Pkt 5 Pairs



Product Dimensions:

4mm Small:	70mm x 56mm approx.
4mm Medium:	76mm x 62mm approx.
4mm Large:	81mm x 64mm approx.
6mm Small:	102mm x 56mm approx.
6mm Medium:	107mm x 62mm approx.
6mm Large:	116mm x 64mm approx.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987