


Hapla

Quality
Adhesive System

Over 35 years of Scientific and Clinical development.





Not all wool felts are the same...

What makes **Hapla** superior?

THE WOOL

Superior selection = Superior performance

- Hapla chooses only the finest grade wool with minimal medulla core structure to ensure ideal pressure re-distribution and absorption behaviour
- Hapla maintain tight control on wool staple length. This ensures consistent bonding of wool fibres, maintains a soft and comfortable felt bed and minimises de-lamination



Acrylate: The building block of the Hapla Adhesive System

THE ADHESIVE

- Function specific formulation applied with the unique Hapla coating technology
- Hypo-allergenic adhesive that is breathable and moisture permeable
- To achieve optimal characteristics the adhesive is energised - allowing chemical cross-bonding with the wool fibres

Images are property and copywrite of CSIRO. www.scienceimage.csiro.au. Wool follicles under the microscope, photographer - Livestock Industries.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com