





HEELUPS™

- Adjustable heel insert to relieve Achilles tendon and calf injuries
- Insert offers a rise of 15mm, 10mm or 5mm allowing the height to be gradually reduced during recovery from injury
- Designed in Australia by leading experienced clinicians using the latest scientific principles
- Includes wedges to attach to the underside of inserts to create a varus/supinated or valgus/pronated position of the heel

HEELUPS-S	Small	👤 3½-7	👤 4½-8	15cm*
HEELUPS-M	Medium	👤 7½-11	👤 8½-13	17cm*
HEELUPS-L	Large	👤 11½-16	👤 ---	18cm*

ARTG DECLARED
407334

*approx. orthotic length. Sizing shown in U.S Sizes



SUPPORTS RECOVERY FROM ACHILLES TENDON & CALF MUSCLE INJURIES

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com