

HEINE DELTAone

The most important points at a glance



The "Application" video shows the DELTAone in action.



The "DELTAone DERM app" video shows you how to use the DELTAone with the exclusive HEINE DERM app for digital documentation purposes.



DELTAone SPECIAL PODIATRY KIT

- Includes autoclavable small contact plate and universal phone connector
- LED^{HQ} technology that evenly distributes illumination from edge-to-edge, without distracting reflections, shadows, bright or dark spots
- True field of view 22mm lens provides 10x magnification
- Switch between polarised and non-polarised illumination mode at the touch of a button
- Fast change from contact to non-contact examination
- 3 light intensity levels
- Pocket-sized for portability
- Extremely sharp image with the new achromatic HEINE optical system
- Free HEINE Derm app available for iOS Phones (model 6 or later)
- Includes rechargeable Lithium-Ion battery with Micro USB cable

DELTAONEKT Heine DELTAone Dermatoscope Podiatry Kit

Log-in to www.briggatemedical.com for pricing



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com