



9	43
7.5	42
6	41
4.5	40
3	39
1.5	38
13	37
11.5	36
10	35
8.5	34
7	33
5.5	32
	31
	30
	29
	28
	27
	26
	25
	24
	23
	22
	21

U.S. KIDS EURO  
**Shoe Size**



**littleSTEPS®**  
foot orthotics for kids  
**FOOT CHART KEY:**  
LS = littleSTEPS® foot orthotics size  
GP = littleSTEPS® gait plate size

PLACE PATIENT'S  
FOOT ON THE LINE

SIZE	MODEL	U.S. KIDS	U.S. MEN	U.S. WOMEN	EURO
9	YOUTH/ADULT*	7.5 to 8.5	7.5 to 8.5	9 to 10	41 to 42
8	YOUTH/ADULT*	6 to 7	6 to 7	7.5 to 8.5	39 to 40
7	YOUTH/ADULT*	4.5 to 5.5	4.5 to 5.5	6 to 7	37 to 38
6	CHILDREN	3 to 4			35 to 36
5	CHILDREN	1.5 to 2.5			33 to 34
4	CHILDREN	13 to 1			31 to 32
3	CHILDREN	11.5 to 12.5			29 to 30
2	CHILDREN	10 to 11			27 to 28
1	TODDLER	8.5 to 9.5			25 to 26
0	TODDLER	7 to 8			23 to 24
00	TODDLER	5 to 6.5			21 to 22

\*Youth sizes are roughly equivalent to Men's sizes in length.

**IMPORTANT: WHEN PRINTING - DO NOT SELECT "FIT TO PAGE" AS THIS WILL EFFECT THE SCALE OF THE CHART AND RESULT IN INCORRECT SIZE FOR PATIENT.**