


MELAG[®]

competency in hygiene



INCLUDES YOUR CHOICE OF DOCUMENTATION OUTPUT

Choose between the MelaFlash PC Storage unit or the traditional MelaPrint Printer.



MAXIMISED CHAMBER CONFIGURATION

Includes a Pouch Rack and Tray to maximise the cycle load



Euroklav[®] 23 S+

22 litre Chamber

- S-Class Steriliser
- Approx. Cycle Time: 45-50 Minutes
- Includes 1 Tray, 1 Pouch Rack and your choice of Data Storage
- Monitors water quality
- Automatic drying cycle
- Stainless Steel chamber for easy cleaning

Specifications:

- External Dimensions:
- 43cmW x 52cmH x 66cmD (22L)
- Weight: 43kg
- Chamber Dimensions:
- 23cmW x 45cmD (22L)
- ARTG# 146443

MELA23S

Log-in to www.briggatemedical.com for pricing

- Melaflash OR Melaprint included
- Supplied in an un-calibrated state
- Calibration/validation is required once the unit reaches it's final destination

Warranty:
2 years / 2,000 cycles

Manufactured in:
Germany



Delivery fees may apply

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com