

# Mepilex®

## Effective absorbent foam dressing

### Polyurethane foam pad

- Absorbs exudate
- Conformable
- May be used under compression

### Polyurethane backing film

- Breathable
- Shower-proof<sup>1</sup>
- Bacteria and viral barrier<sup>2,3,†</sup>

†Microbes larger than 25 nm



Safetac  
TECHNOLOGY



### Safetac® wound contact layer

- Less pain during dressing changes<sup>4</sup>
- Does not stick to the wound, for less disturbance<sup>5</sup>
- Seals wound margins to avoid maceration<sup>6</sup>



## How Mepilex works

Mepilex is a soft and highly conformable foam dressing that absorbs exudate and maintains a moist wound environment.

The Safetac technology layer seals the wound edges, preventing the exudate from leaking onto the surrounding skin, thus minimising the risk of maceration. The Safetac technology layer ensures that the dressing can be changed without damaging the wound or surrounding skin or exposing the patient to additional pain.

## Benefits of Mepilex

- Minimises pain and trauma at dressing changes<sup>1,4,5</sup>
- Stays in place allowing easy application of compression or retention bandages
- Well-suited for use under compression bandages
- May be cut to suit various wound shapes and difficult-to-dress locations
- Promotes patient comfort during wear
- May remain in place for several days depending on the condition of the wound
- May be adjusted without losing its adherent properties
- Low potential for skin irritation and allergy

## Areas of use

You can use Mepilex to treat a wide range of exuding wounds such as:

- Moderately exuding wounds
- Leg and foot ulcers
- Pressure injuries
- Under compression
- Graft and donor sites
- Traumatic wounds (e.g. skin tears and wounds healing by secondary intention)

## Mepilex ordering information†

Product Code	Size	Pcs/box
294015	5 x 5 cm	5
294100	10 x 10 cm	5