

Mepilex® XT

Absorbent foam dressing with exudate channels

Polyurethane foam pad

- Enhanced exudate management capability
- Conformable and flexible foam

Breathable polyurethane backing film

- Shower-proof[†]
 - Barrier to virus and bacteria^{2,3,†}
- † Microbes larger than 25 nm



Safetac®
TECHNOLOGY

Safetac® technology layer

- Less pain during dressing changes⁴
- Does not stick to the wound, for less disturbance⁵
- Seals wound margins to avoid maceration⁶



Handles both low and high viscosity fluid^{7,8}

Ultra-absorbent foam with integrated channels absorbs exudate, drawing it away from the wound bed, then spreading and absorbing it across the foam pad.



How Mepilex XT works

Mepilex XT is a soft and conformable⁹ foam dressing with integrated channels, designed to transfer exudate rapidly⁸ into the absorbent foam pad. Safetac technology seals the wound edges, preventing exudate from leaking onto the surrounding skin, thus minimising the risk of maceration. Safetac also ensures that the dressing can be changed with minimised risk of damaging the wound or surrounding skin or exposing the patient to additional pain^{10,11}.


Benefits of Mepilex XT

- Minimises maceration, absorbs both low and high viscosity exudate^{7,8}
- Can be used in all exuding wound healing stages
- Minimises pain and trauma at dressing changes^{4,5}
- Promotes patient comfort
- Well-suited for use under compression bandages^{7,8}
- Can remain in place for several days depending on the condition of the wound^{9,12}
- Can be cut to suit various wound shapes and difficult-to-dress locations
- Stays in place allowing for "hands-free" application of compression or retention bandages⁸
- Can be lifted and adjusted without losing its adherent properties
- Non-sensitising¹³


Areas of use

You can use Mepilex XT to treat a wide range of exuding wounds such as:

- Exuding acute and chronic wounds in all healing stages
- Leg and foot ulcers
- Pressure injuries
- Traumatic wounds



Pressure injury



Foot ulcer

Mepilex XT ordering information[‡]

Product Code	Size	Pcs/box
211015	5 x 5 cm	5
211100	10 x 10 cm	5