

Milpara Podiatry Cabinet



CABINETWRT
White with Stainless Steel Top



CABINETSR
Silver with Stainless Steel Top

NOW
INCLUDES
SOFT-CLOSE
DRAWERS

Milpara Podiatry Cabinet

The Milpara Podiatry Cabinet has castors, two large drawers, two small drawers, and a glide out tray to suit most podiatry drills. All drawers are on runners and are fully melamine lined for easy cleaning. Supplied fully assembled and ready to use.


CABINETWR	White with Drill Tray on Right
CABINETWL	White with Drill Tray on Left
CABINETSR	Silver with Drill Tray on Right
CABINETSL	Silver with Drill Tray on Left
CABINETWRT	White with Drill Tray on Right & Stainless Steel Top
CABINETWLT	White with Drill Tray on Left & Stainless Steel Top
CABINETSR	Silver with Drill Tray on Right & Stainless Steel Top
CABINETSLT	Silver with Drill Tray on Left & Stainless Steel Top

Log-in to www.briggatemedical.com for pricing



- Australian laminate
- Australian aluminum handle profile
- Australian stainless steel
- Australian castors
- European SOFT-CLOSE sliders
- Handcrafted by Peter (Melbourne)

Delivery fees may apply depending on delivery destination and courier service required.



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com