

Facehole Headrest



Multi Purpose Patient Chair

- Electric Lift and Seat Tilt functions via foot controls
- Adjustable gas-strut Legrest & Backrest sections
- Flip-up armrests
- Face-hole with headrest plug
- Debris tray
- Retractable castors
- Adjustable feet for uneven surfaces
- Heavy duty PVC vinyl upholstery

250 KG
LIFT CAPACITY

WARRANTY:
3
YEARS

FREIGHT
CHARGES
MAY APPLY



CHAIR14-NB Navy Blue

SPECIFICATIONS:	
Max. Lift Capacity	250kgs
Minimum Seat Height	56.5cm
Maximum Seat Height	98cm
Motorised seat tilt	+25°
Maximum Length (flat)	220cm
Width (including armrests)	91cm
Seat Cushion	60cmW x 56cmD
Back Cushion	60cmW x 91cmH
Leg Cushion	60cmW x 38cmD
Cushion Thickness	5cm
Cushion Density	36kg/m3

◆ PRODUCT PRICING CAN BE FOUND ON OUR WEBSITE

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com