

NAMROL Grindmaster PN-203

Warranty:
1 year

Delivery fees
may apply



Namrol Grindmaster PN-203

The Namrol PN-203 GrindMaster is ideal for grinding orthotics made of leather, foam, plastic or cork, as well as for finish polishing of plastic materials. It can be easily positioned on a counter, bench, or table-top to fit almost anywhere in your laboratory.

Features:

- Dual Wheel Grinder
- Fully integrated bi-directional dust extraction system
- Full 3/4Hp, 2950 RPM
- Chuck Mechanism permits easy changing of sanding drums
- Replaceable dust bag for easy, convenient cleaning

Includes:

- 1 x fitted Dust Collection Bag & 2 x spare Dust Bags
- 1 x 10cm drum with 36 Grit Coarse Sanding Belt
- 1 x 80 Grit Fine Sanding Belt
- 1 x 30mm Finishing Wheel
- 1 x 30mm Expanding Drum

Specifications:

- Dimensions: 450mmH x 700mmL x 630mmD
- Weight: 37kg

NAMROL Namrol Grindmaster PN-203
Benchtop Clinic Grinder with Dust Extraction.

Log-in to www.briggatemedical.com for pricing



Additional accessories are available for
the Namrol GrindMaster.
Contact **Briggate Medical** for
further information.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggatemedical.com

www.briggatemedical.com