

Operators Chairs



Nickel Grey



Royal Blue



Anthracite



Ocean Teal



Lapis Blue



Platin Grey



Fuchsia

'V' Series Operators Chair

- Hydraulic ram height adjustment
- Polished chrome cast aluminium star base
- Friction lock castors
- Double layer Dunlop foam seat cushion with anti-microbial treatment

- Minimum Height: 450mm
- Maximum Height: 600mm
- Cushion Dimensions: 420mmD x 70mmH
- Safe Working Load: 110kg

Log-in to www.briggatemedical.com for pricing

- OPCHAIRV01 Platin Grey
- OPCHAIRV02 Royal Blue
- OPCHAIRV04 Anthracite
- OPCHAIRV06 Lapis Blue
- OPCHAIRV08 Black
- OPCHAIRV09 Navy
- OPCHAIRV10 Midnight Blue
- OPCHAIRV11 Slate Grey
- OPCHAIRV12 Ocean Teal
- OPCHAIRV13 Fuchsia



on Platin Grey, Royal Blue, Anthracite, Ocean Teal & Fuchsia chairs

Colour matched to suit most Patient Chairs available from Briggate Medical.

Delivery fees may apply.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987