

# Ortho-Press



## Ortho-Press Vacuum Press

- 4-5 second pull down
- Easy to use controls
- Vacuum hold control
- Hassle free diaphragm replacement

### Dimensions:

- 230mmH x 500mmL x 310mmW
- Weight: 18kgs

**PRESS** Ortho-Press Vacuum Press

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing

Warranty:  
**1 year**

Manufactured in:  
**Australia**



Delivery fees  
may apply

## Ortho-Press "SLAVE" Vacuum Press

The Ortho-Press "Slave" Vacuum Press is available as an optional extra. It connects directly to the main Ortho-Press Vacuum Press using the pump from that unit so that two orthotics can be formed simultaneously.

**PRESSLAV** Ortho-Press "SLAVE" Vacuum Press **Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing**

### Replacement Diaphragm

<b>DIAPH01</b>	Pre-Cut	400mm x 560mm	1.59mm
<b>DIAPH02</b>	Blank	450mm x 610mm	1.59mm



# briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

[briggate@briggate.com.au](mailto:briggate@briggate.com.au)

[www.briggatemedical.com](http://www.briggatemedical.com)