


PPT403 Heel Wedges

SPECIFICATION SHEET



PPT403 Heel Wedges SMALL

Packet of 6 pairs
Manufacturers Code: 22200
Self-adhesive backing

Approx. Width: 63mm
Approx. Length: 60mm
Approx. Peak Height: 5mm

PPT403 Heel Wedges MEDIUM

Packet of 6 pairs
Manufacturers Code: 22210
Self-adhesive backing

Approx. Width: 69mm
Approx. Length: 67mm
Approx. Peak Height: 5mm

PPT403 Heel Wedges LARGE

Packet of 6 pairs
Manufacturers Code: 22220
Self-adhesive backing

Approx. Width: 79mm
Approx. Length: 79mm
Approx. Peak Height: 5mm

Constructed from PPT open-cell polyurethane foam.

Key properties:

- retains it's shape
- provides comfort and shock absorption
- breathable to keep skin cool and dry

Product outlines shown above are approximate to scale.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987