


PPT404 Heel Lifts SPECIFICATION SHEET



PPT404 Heel Lifts SMALL

Packet of 6 pairs
Manufacturers Code: 22300
Self-adhesive backing
Approx. Width: 60mm
Approx. Length: 63mm
Approx. Peak Height: 6.5mm

PPT404 Heel Lifts MEDIUM

Packet of 6 pairs
Manufacturers Code: 22310
Self-adhesive backing
Approx. Width: 66mm
Approx. Length: 65mm
Approx. Peak Height: 6.5mm

PPT404 Heel Lifts LARGE

Packet of 6 pairs
Manufacturers Code: 22320
Self-adhesive backing
Approx. Width: 74mm
Approx. Length: 80mm
Approx. Peak Height: 6.5mm

Constructed from PPT open-cell polyurethane foam.
Key properties:
- retains it's shape
- provides comfort and shock absorption
- breathable to keep skin cool and dry

Product outlines shown above are approximate to scale.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987