


## PPT406 Metatarsal Bars SPECIFICATION SHEET



### PPT406 Metatarsal Bars SMALL

Packet of 6 pairs  
Manufacturers Code: 22500  
Self-adhesive backing  
Approx. Width: 60mm  
Approx. Length: 60mm  
Approx. Peak Height: 5mm

### PPT406 Metatarsal Bars MEDIUM

Packet of 6 pairs  
Manufacturers Code: 22510  
Self-adhesive backing  
Approx. Width: 67mm  
Approx. Length: 65mm  
Approx. Peak Height: 5mm

### PPT406 Metatarsal Bars LARGE

Packet of 6 pairs  
Manufacturers Code: 22520  
Self-adhesive backing  
Approx. Width: 76mm  
Approx. Length: 69mm  
Approx. Peak Height: 5mm

Constructed from PPT open-cell polyurethane foam.  
Key properties:  
- retains it's shape  
- provides comfort and shock absorption  
- breathable to keep skin cool and dry

Product outlines shown above are approximate to scale.

# briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987