


PPT407 Metatarsal Pads SPECIFICATION SHEET



PPT407 Metatarsal Pads SMALL

Packet of 6 pairs
Manufacturers Code: 22600
Self-adhesive backing
Approx. Width: 53mm
Approx. Length: 48mm
Approx. Peak Height: 6mm

PPT407 Metatarsal Pads MEDIUM

Packet of 6 pairs
Manufacturers Code: 22610
Self-adhesive backing
Approx. Width: 54mm
Approx. Length: 65mm
Approx. Peak Height: 7mm

PPT407 Metatarsal Pads LARGE

Packet of 6 pairs
Manufacturers Code: 22620
Self-adhesive backing
Approx. Width: 60mm
Approx. Length: 65mm
Approx. Peak Height: 8mm

Constructed from PPT open-cell polyurethane foam.
Key properties:
- retains it's shape
- provides comfort and shock absorption
- breathable to keep skin cool and dry

Product outlines shown above are approximate to scale.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987