


PPT462 Heel Spur Inserts

SPECIFICATION SHEET



PPT462 Heel Spur Insert SMALL

Packet of 6 pairs
Manufacturers Code: 13801
Self-adhesive backing

Approx. Width: 55mm
Approx. Length: 70mm
Approx. Peak Height: 3mm

PPT462 Heel Spur Insert LARGE

Packet of 6 pairs
Manufacturers Code: 13803
Self-adhesive backing

Approx. Width: 70mm
Approx. Length: 80mm
Approx. Peak Height: 3mm

Constructed from PPT open-cell polyurethane foam.
Key properties:

- retains it's shape
- provides comfort and shock absorption
- breathable to keep skin cool and dry

Product outlines shown above are approximate to scale.

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987