

Peg-Assist™ Insole

Off-loading Insole System



The DARGO Peg-Assist™ System features a multi-indication removable-peg chasis that effectively off-loads the plantar aspect of the foot after surgery or when wounds are present.


DARGO




Peg-Assist™ Features and Benefits:

- › Removable pegs allow for localized off-loading of wounds and ulcerations of the foot
- › Pressure reduced by as much as 60% while allowing patient to remain ambulatory.
- › Poron® cover eliminates the incidence of ring edema and edge abrasion.
- › Included stabilizer board prevents adjacent pegs from collapsing.
- › Available in a wide range of sizes and configurations designed to fit most DARCO shoes and walkers


Peg-Assist Insole™ Cross Section



Focal pressure reduction of 60% under the 1st MPJ is achieved by removing the pegs beneath the area marked by the red circle.



Standard Surgical Shoe



DARCO Med-Surg™ Shoe with Peg-Assist™ Insole

Insole for Square Toe Med-Surg Shoe (7 sizes)

Item Number Women's	Item Number Men's	Description	Size
PTQW1	PTQM1	Peg-Assist	Small
PTQW2	PTQM2	Peg-Assist	Medium
PTQW3	PTQM3	Peg-Assist	Large
-	PTQM4	Peg-Assist	X-Large

Insole for Other DARCO Square Toe Shoes and Walkers (5 sizes)

Item Number	Description	Size
PQ0	Peg-Assist	X-Small
PQ1	Peg-Assist	Small
PQ2	Peg-Assist	Medium
PQ3	Peg-Assist	Large
PQ4	Peg-Assist	X-Large

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA

Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037

info@darcointernational.com | www.darcointernational.com

As with all DARCO products, a portion of the profits from the sale of these items will be donated to encourage and support continuing education in our profession.