

Med-Surg™ Shoe

Post Operative Protective Footwear



The DARCO Med-Surg™ Shoe is the most versatile post operative shoe on the market today. It's fully customizable insole is double-padded and removable to accommodate the Peg-Assist™ off-loading insole.

DARCO



Med-Surg™ Features and Benefits:

- › **Square Toe Design** acts as a bumper to provide extra protection and facilitates better left/right fit.
- › **MetaShank™** protection provides more rigid control under the metatarsal heads.
- › **Stapless Closure** eliminates buckle pressure on the forefoot and can expand to accommodate dressing up to 17 inches in circumference.
- › **Rocker Sole** reduces plantar pressure on the heel and forefoot by **over 25%**.
- › Fully compatible with **DARCO's Peg-Assist™ off-loading insole**.

Med-Surg™ Indications:

- › After foot surgery or when K-wires are present
- › When large bandages or dressing are required
- › As a platform for the Peg-Assist™ insole when off-loading wounds

Women's Item Number	Men's Item Number	Description	Size	Color	Women's Size	Men's Size
MQW1B	MQM1B	Med-Surg	X-Small	Black	4 - 6	6 - 8
MQW2B	MQM2B	Med-Surg	Small	Black	6.5 - 8	8.5 - 10
MQW3B	MQM3B	Med-Surg	Medium	Black	8.5 - 10	10.5 - 12
	MQM4B	Med-Surg	Large	Black		12.5 - 14

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA

Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037

info@darcointernational.com | www.darcointernational.com

As with all DARCO products, a portion of the profits from the sale of these items will be donated to encourage and support continuing education in our profession.