

PowerKids®

Affordable, easy-to-use, prefabricated orthotics for kids – similar fit and function of custom-made products. Provides rear-foot stability by properly positioning and cradling the heel. Built-in arch support for kids vertically lifts and supports the arch and restores proper alignment to the mid-foot and forefoot. These 3/4 length kids insoles are comfortable on the forefoot and toes.

Materials:

- Fabric top cover with anti-microbial protection, Royal Blue with Paediatric Print
- Poron® foam top layer, Blue, 1.6mm approx.
- EVA base, Red

Features and Benefits:

- Built-in arch support
- Integrated 2° medial heel post and deep heel cup for hind-foot stability and motion control
- 3/4 length can be easily transferred from shoe-to-shoe
- Cushioning, shock-absorbing materials
- Friction reducing, anti-microbial top fabric

Indications:

- Mild to moderate over-pronation (flat feet), lacking arch and foot support, lacking proper alignment and balance, excessive motion and muscle fatigue, foot pain and discomfort

POWERKID-B	US: 8½-9½	10.5cm*
POWERKID-C	US: 10-11	11cm*
POWERKID-D	US: 11½-12½	12cm*
POWERKID-E	US: 13-1	12.5cm*
POWERKID-F	US: 1½-2½	13cm*
POWERKID-G	US: 3-4	14cm*

*approx. orthotic length.



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com