

Powerstep

PROTECH - CONTROL

Firmer arch support, maximum cushioning and a 2-degree heel wedge for greater stability and correction for overpronation or fallen arches.



POW-01AS	Size AS	3-3½	5-5½	23.5cm*
POW-01A	Size A	4-4½	6-6½	24.5cm*
POW-01B	Size B	5-5½	7-7½	25cm*
POW-01C	Size C	6-6½	8-8½	25.5cm*
POW-01D	Size D	7-7½	9-9½	26cm*
POW-01E	Size E	8-8½	10-10½	27cm*
POW-01F	Size F	9-9½	11-11½	28cm*
POW-01G	Size G	10-10½	12-12½	29cm*
POW-01H	Size H	11-11½	12-12½	29.5cm*
POW-01I/J	Size I/J	12-13	---	30cm*
POW-01K	Size K	14-15	---	32cm*

*approx. orthotic length. U.S Sizing shown.

- Strong, medical-grade foot support calibrated with flexibility
- Deeper heel cup and 2-degree medial post for improved foot control and positioning
- Stronger arch firmness
- Double layer, shock absorbing EVA/VCT® cushion casing for enhanced comfort
- Friction reducing anti-microbial top fabric
- Poron® heel pad cushions the heel

INDICATIONS:

Moderate Over-Pronation; Plantar Fasciitis/Heel Pain; Foot, Arch and Heel Pain; Calluses; Hammertoes and Claw Toe; Bunions; Neuromas; Ankle and Knee Pain

IDEAL SHOE TYPE:

Athletic Shoes; Roomier Dress and Casual Shoes

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com