


Powerstep

PROTECH - CONTROL 3/4

The ProTech Control 3/4 length shoe orthotics are versatile enough to fit in a variety of shoes. The firmer but flexible arch support and deeper heel cup with 2-degree medial post provide improved foot control and positioning.



POW-03AS	Size AS	4-4½	6-6½	16cm*
POW-03A	Size A	5-6½	7-8½	17cm*
POW-03B	Size B	7-8½	9-10½	18cm*
POW-03C	Size C	9-9½	11-11½	18.5cm*
POW-03D	Size D	10-10½	12-12½	19cm*
POW-03E	Size E	11-11½	13-13½	19.5cm*
POW-03F	Size F	12-13½	14-15½	20cm*

*approx. orthotic length. U.S Sizing shown.

- Strong, medical-grade foot support calibrated with flexibility
- Deeper heel cup and 2 degree medial post for improved foot control and positioning
- Stronger arch firmness
- Double layer, shock absorbing EVA/VCT® cushion casing for enhanced comfort
- Friction reducing anti-microbial top fabric
- Poron® heel pad cushions the heel

INDICATIONS:

Over-Pronation (Mild to Moderate); Plantar Fasciitis/Heel Pain; Foot, Arch and Heel Pain; Calluses; Hammertoes and Claw Toe; Bunions; Neuromas; Ankle and Knee Pain

IDEAL SHOE TYPE:

Dress and Casual Shoes, Slip-On Dress Shoes, Low to Mid Heel Dress Shoes, Athletic Cleats

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com