


Powerstep

WIDE

PROTECH - CONTROL

Maximum cushioning, firmer arch support, and 2-degree heel wedge for greater support and to prevent "bottoming out" under weight. Wider toebox fits shoes 3E-6E.



POW-02A	Size A	4-4½	6-6½	24.5cm*
POW-02B	Size B	5-5½	7-7½	25cm*
POW-02C	Size C	6-6½	8-8½	25.5cm*
POW-02D	Size D	7-7½	9-9½	26cm*
POW-02E	Size E	8-8½	10-10½	26.5cm*
POW-02F	Size F	9-9½	11-11½	27.5cm*
POW-02G	Size G	10-10½	12-12½	29cm*
POW-02H	Size H	11-11½	12-12½	29.5cm*
POW-02I/J	Size I/J	12-13	---	30cm*
POW-02K	Size K	14-15	---	32cm*

*approx. orthotic length. U.S Sizing shown.

- Semi-rigid support shell calibrated for firmer arch support to resist "bottoming out" under weight
- Deeper heel cup and 2-degree medial post for improved foot control and positioning
- Wider profile to better fit wide-width shoes, 3E - 6E
- Double-layer, low-profile cushioning for enhanced comfort
- Friction reducing anti-microbial top fabrica

INDICATIONS:

Those Needing Firmer, Heavier-Weight Support; Wide Feet; Moderate Pronation; Plantar Fasciitis/Heel Pain; Foot, Arch and Heel Pain; Calluses; Hammertoes and Claw Toe; Bunions; Ankle and Knee Pain; Metatarsalgia

IDEAL SHOE TYPE:

Wide Shoes from 3E-6E, Comfort or Athletic Shoes, Therapeutic Shoes, Casual Lace Up Shoes

briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com