



# ri-san<sup>®</sup> aneroid sphygmomanometer

The ri-san high-quality, palm-style sphygmomanometer is the embodiment of ergonomic design. The air release valve is conveniently regulated by finger-tip control.

- Specially tempered copper-beryllium diaphragm, almost non-aging
- Precision, wear-free air release valve offers fine adjustment
- Microfilter protects air release valve and measuring system
- Diaphragm can be pressurized up to 600 mmHg
- Latex-free Cuff and Bulb
- Tube connection at the top of the pressure gage for unhindered measurement
- Suitable for right or left handed users
- Easily legible scale: Ø 64 from 0 to 300 mmHg
- Maximum error tolerance of +/- 3 mmHg

## Riester Ri-san Sphygmomanometer

Package includes:

- 1 x Riester Ri-san Sphygmomanometer
- 1 x Adult Cuff (size range: 24cm-34cm)
- 1 x Vinyl zip-up protective case

**RI-SAN**

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing



## Riester Ri-san Sphygmomanometer with Hokanson 2.5cm Digital Cuff

Package includes:

- 1 x Riester Ri-san Sphygmomanometer
- 1 x Adult Cuff (size range: 24cm-34cm)
- 1 x Digital Cuff (2.5cm)
- 1 x Vinyl zip-up protective case

**RI-SAN-TBI**

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing



INCLUDES DIGITAL CUFF FOR CONDUCTING TOE BRACHIAL INDEX MEASUREMENTS\*

\*in conjunction with a doppler that has T.B.I. capabilities.

## Riester Ri-san Sphygmomanometer Multi-Cuff Kit

Package includes:

- 1 x Riester Ri-san Sphygmomanometer
- 1 x Infant Cuff (size range: 8-13cm)
- 1 x Small Adult Cuff (size range: 17-24cm)
- 1 x Adult Cuff (size range: 24cm-34cm)
- 1 x XXL Adult Cuff (size range: 32-41cm)
- 1 x Vinyl zip-up protective case

**RI-BPMCKIT**

Log-in to [www.briggatemedical.com](http://www.briggatemedical.com) for pricing



# briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

[briggate@briggate.com.au](mailto:briggate@briggate.com.au)

[www.briggatemedical.com](http://www.briggatemedical.com)