

SaniGrinder704

Warranty:
1 year

Manufactured in:
USA

Delivery fees
may apply



Sani Grinder 704 Clinic Grinder

- Fully integrated dust extraction system
- 2 Speed 1/4hp Red Wing Motor
- Chuck Remover permits easy changing of sanding drums
- Suction Hood rotates to suit your grinding and polishing needs
- Includes: 1 x Dust Collection Bag, 1 x 80 Grit Medium Sanding Belt

Dimensions:


- 250mmH x 630mmL x 340mmD
- Drum: 100mm x 80mm
- Weight: 20.5 kg

SANI-704

Sani Grinder 704

Benchtop Clinic Grinder with Dust Extraction.

Log-in to www.briggatemedical.com for pricing



An extensive range of Sanding Belts, Drums, Polishing buffs and other accessories are also available from Briggate Medical.



briggatemedicalcompany

Supporting Podiatrists with integrity and trust since 1987

(03) 8586 7800 (Melb. Metro)

1800 33 4142 (Free Call)

briggate@briggate.com.au

www.briggatemedical.com